

AIRWRAP™ COMPRESSION BANDAGE

Hands-free direct pressure



Keeps direct pressure on a wound from the point of injury to definitive care

- Provides focal, adjustable pressure
- 50% reduction in treatment time and 75% reduction in blood loss compared to traditional methods*
- Ideal for tourniquet conversions

Using AIRWRAP, responders are freed to treat additional injuries or patients

- Ideal for multi-site, multi-patient triage in mass-casualty events

Capable of treating a wide variety of bleeding injuries almost anywhere on the body

Standard issue in Medic/Corpsmen Kits

*Kragh, J. F. et al. (2014). A manikin model for study of wound-packing interventions to control out-of-hospital hemorrhage. The American Journal of Emergency Medicine, 32(9), 1130-1131.

PRODUCT SPECIFICATIONS

AIRWRAP
NSN 6530-01-643-5297



KEY ATTRIBUTES

4"compression bandage + bladder
Pop-up inflation indicator
Can be filled with IV fluid using luer lock or 60cc syringe
5 year shelf life

PRODUCT DIMENSIONS

LENGTH: 152.4 cm (60 in)
WIDTH: 10 cm (4 in)

PACKAGED DIMENSIONS

9.5 x 12.7 x 5.1 cm
(3.7 x 5 x 2 in)
WT: 0.14 kg (0.3 lbs)

AIRWRAP XL
NSN 6530-01-657-1978



KEY ATTRIBUTES

6"compression bandage + removable bladder (for additional treatment options)
Can be filled with IV fluid using luer lock or 60cc syringe
5 year shelf life

PRODUCT DIMENSIONS

LENGTH: 219.7cm (86.5 in)
WIDTH: 15 cm (6 in)

PACKAGED DIMENSIONS

10.75 x 15.24 x 5.1 cm
(4.2 x 6 x 2 in)
WT: 0.23 kg (0.5 lbs)