

PARABELT® TOURNIQUET BELT

Be prepared anytime anywhere to stop bleeding and save a life



Daily-wear tourniquet belt combination designed for immediate hemorrhage control in any situation

Lightweight, compact and comfortable

Simple to use



- Intuitive ratcheting design enables fast application with minimal training
- Thumb loop doubles as bite strap for one-handed application

Safe & Effective

- Uses same CoTCCC-RECOMMENDED tightening system as TX Tourniquets
- Precise, one-handed micro-adjustments to treat re-bleeding and minimize patient discomfort

Ideal for military and civilian first responders, law enforcement, outdoors enthusiasts or concerned citizens

PRODUCT SPECIFICATIONS

PARABELT MODELS	KEY ATTRIBUTES	PACKAGED DIMENSIONS
BLACK 	Lightweight Type 13 MILSPEC webbing belt with integrated 1.75 in wide tourniquet Safety rated at 4000 lb	WT: 0.54 kg (2.02 lbs)
TAN 	Mechanical ratchet with protected release Rated 18kN AustriaAlpin COBRA® frame buckle Best performance in study of belts used as tourniquets*	
SIZE CHART		
Size	Waist (in)	
S	26-30	
M	30-35	
L	35-40	
XL	39-44	
XXL	44-50	

* Bequette, B. W., Kragh, J. F., Aden, J. K., & Dubick, M. A. (2017). Belts Evaluated as Limb Tourniquets: BELT Study Comparing Trousers Used as Medical Devices in a Manikin Model of Wound Bleeding. *Wilderness & Environmental Medicine*, 28(2), 84-93.