

TX™ TOURNIQUETS

Intuitive to use and made to last



CoTCCC RECOMMENDED

Available in 2" (TX2) and 3" (TX3) widths

- TX3 is the widest non-pneumatic tourniquet on CoTCCC-recommended list

Intuitive ratcheting design enables fast application with minimal training

- Self-locking, non-slip ratchet
- Loop assist enables one-handed application

Durable construction allows repeated use in training without loss of function

- No need to purchase additional training tourniquets - train with the same TQ you take to combat
- High-density polymer ladder maintains function in presence of sand and dirt
- UV-resistant MILSPEC webbing
- No expiration date

TX design allows precise, micro-adjustments to treat rebleeding and minimize patient discomfort

PRODUCT SPECIFICATIONS

TX2
NSN 6515-01-667-6027



KEY ATTRIBUTES

2 in (5 cm) webbing width
Ratchet design with mechanical
release

PRODUCT DIMENSIONS

LENGTH: 152.4 cm (60 in)
WIDTH: 5.1 cm (2 in)

PACKAGED DIMENSIONS

21.6 x 5.1 x 3.8 cm
(8.5 x 2 x 1.5 in)

WT: 0.09 kg (0.2 lbs)

TX3
NSN 6515-01-667-6208



KEY ATTRIBUTES

3 in (7.6 cm) webbing width
Ratchet design with mechanical
release
Designed for Prolonged Field Care
where tissue and nerve damage
are a major concern

PRODUCT DIMENSIONS

LENGTH: 152.4 cm (60 in)
WIDTH: 7.6 cm (3 in)

PACKAGED DIMENSIONS

21.6 x 7.6 x 3.8 cm
(8.5 x 3 x 1.5 in)

WT: 0.18 kg (0.4 lbs)