

XGAUZE® WOUND DRESSING

Expands on contact to simplify wound packing and control hemorrhage



The only gauze dressing that absorbs, expands and applies pressure

- Equivalent to 3 rolls of standard gauze in absorption and expansion
- Utilizes similar compressed sponge technology as CoTCCC recommended XSTAT
- Applies over 50% higher pressure within a wound than standard gauze*

Simplifies wound packing

- Alleviates need for wound-packing expertise
- Does not require manual compression, unlike hemostatic gauze
- Reduces application time by over 40% vs standard gauze*

Standard issue in Medic/Corpsmen Kits

Hundreds of successful, documented combat uses

*Bullard, M. et al. (2016, August). *Comparison of XGAUZE® Trauma Dressing to Standard Gauze in a Benchtop Wound Model*. Poster session presented at Military Health System Research Symposium, Kissimmee, FL.

PRODUCT SPECIFICATIONS

XGAUZE
NSN 6510-01-643-5391



KEY ATTRIBUTES

Z-folded, sterile gauze embedded with
~138 expanding minisponges

Compressed sponge technology
similar to XSTAT

Radiopaque marker sewn into dressing

Absorbs up to 400 mL fluid

Vacuum-sealed in printed package

5 year shelf life

PRODUCT DIMENSIONS

LENGTH: 91.44 cm (36 in)

WIDTH: 7.5 cm (3 in)

PACKAGED DIMENSIONS

8.9 x 12.7 x 2.5 cm
(3.5 x 5 x 0.98 in)

CUBE: 0.29 L

WT: 0.05 kg (0.11 lbs)