


XSTAT[™] HEMOSTATIC DEVICE

Revolutionizing the treatment of gunshot and knife wounds



Indicated for treatment of junctional wounds in the groin or axilla and narrow entrance extremity wounds in the arms or legs, for adults and adolescents

CoTCCC RECOMMENDED

The only FDA-cleared device indicated for severe bleeding from penetrating, narrow-entrance wounds

- Applies hemostatic pressure inside the wound cavity to stop bleeding at its source
- Overcomes limitations of traditional wound dressings in treating deep, narrow-tract wounds

Applies in seconds and stops bleeding immediately without manual compression

- Remove from packaging and apply to wound in 30 seconds. XSTAT sponges immediately expand, apply 360° hemostatic pressure and stop bleeding
- No need for the minutes of manual compression necessary with hemostatic gauze

Syringe applicator is easy to learn and easy to use.

- Requires minimal training
- Alleviates need for wound-packing expertise


Compact & lightweight: minimal cube space required

Standard issue in:

SOCOM IFAKs
CASEVAC Sets
Medic/Corpsmen Kits

PRODUCT SPECIFICATIONS

XSTAT 30
Gen 2



Single Pack
NSN 6510-01-644-7335

KEY ATTRIBUTES

CoTCCC Recommended

For larger exit wounds from gunshots or other penetrating trauma

~108 radiopaque sponges per applicator

Vacuum sealed; sterile

5 year shelf life

PRODUCT DIMENSIONS

LENGTH: 14.1 cm (5.55 in)

DIA: 3.0 cm (1.18 in)

PACKAGED DIMENSIONS

17.5 x 28.2 x 3.3 cm

(6.9 x 11.1 x 1.3 in)

CUBE: 0.34 L

WT: 0.10 kg (0.23 lbs)

XSTAT 12



Single Pack
NSN 6510-01-665-5192

3-Pack
NSN 6510-01-657-4737

KEY ATTRIBUTES

CoTCCC Recommended

For smaller entrance wounds from stabbings, shrapnel or small caliber weapons

~38 radiopaque sponges per applicator

Vacuum sealed; sterile

5 year shelf life

PRODUCT DIMENSIONS

LENGTH: 17.5 cm (6.89 in)

DIA: 1.2 cm (0.47 in)

PACKAGED DIMENSIONS

17.5 x 28.2 x 3.3 cm

(6.9 x 11.1 x 1.3 in)

1 Pack

CUBE: 0.32 L

WT: 0.06 kg (0.13 lbs)

3 Pack

CUBE: 0.80 L

WT: 0.13 kg (0.28 lbs)