

Teaches proper use of XSTAT and other hemostatic devices to reduce blood loss in the field



Multi-task training manikin

- Can be used dry or with fluid to simulate bleeding scenarios
- Used to demonstrate effective wound-packing skills and tourniquet/hemostatic device application
- Ideal for crawl, walk, run training
- Designed with three realistic penetrating injuries and sliding tube inserts for easy removal of XSTAT training sponges from wounds

XSTAT 12 and XSTAT 30 Gen 2 training applicators

- Refillable, multi-use applicators with blue tips

XSTAT training sponges

- Blue, non-sterile training sponges are designed for use with XSTAT training manikins and other simulated hemorrhage control training devices

PRODUCT SPECIFICATIONS

Multi-Task Training Manikin NSN 6910-01-668-2607

Two junctional wounds: 1 exit / 1 entry

One mid-thigh wound for gauze packing or tourniquet application

Sliding tube inserts for easy removal of XSTAT sponges from wound

Hoses fitted with luer lock attachment to accept 60cc luer lock syringe or IV line


Detachable hose connections compatible with simulated fluid refill tanks



XSTAT Training Kit NSN 6910-01-668-1798

Kit contains:

- (2) XSTAT 30 training applicators
- (3) XSTAT 12 training applicators
- (1) XSTAT 30 loading funnel
- (1) XSTAT 12 loading tool
- (1) bag of blue training sponges (1000)
- (1) Instruction for Use document



XSTAT Training Sponge Refill Kit NSN 6910-01-667-9079

Kit contains:

- 4000 training sponges in four 1000-count individual bags



WARNING: XSTAT training sponges are non-sterile, do not contain a radiopaque marker and are not for human use.