

# TX™ TOURNIQUETS

Intuitive to use and made to last



## CoTCCC RECOMMENDED

**Available in 2" (TX2) and 3" (TX3) widths**

- TX3 is the widest non-pneumatic tourniquet on CoTCCC-recommended list

**Intuitive ratcheting design enables fast application with minimal training**

- Self-locking, non-slip ratchet
- Loop assist enables one-handed application

**Durable construction allows repeated use in training without loss of function**

- No need to purchase additional training tourniquets - train with the same TQ you take to combat
- High-density polymer ladder maintains function in presence of sand and dirt
- UV-resistant MILSPEC webbing
- No expiration date

**TX design allows precise, micro-adjustments to treat rebleeding and minimize patient discomfort**

## PRODUCT SPECIFICATIONS

**TX2**  
NSN 6515-01-667-6027



### KEY ATTRIBUTES

2 in (5 cm) webbing width  
Ratchet design with mechanical  
release

### PRODUCT DIMENSIONS

LENGTH: 99.1 cm (39 in)  
WIDTH: 5.1 cm (2 in)

### PACKAGED DIMENSIONS

21.6 x 5.1 x 3.8 cm  
(8.5 x 2 x 1.5 in)

WT: 0.09 kg (0.2 lbs)

**TX3**  
NSN 6515-01-667-6208



### KEY ATTRIBUTES

3 in (7.6 cm) webbing width  
Ratchet design with mechanical  
release  
Designed for Prolonged Field Care  
where tissue and nerve damage  
are a major concern

### PRODUCT DIMENSIONS

LENGTH: 152.4 cm (60 in)  
WIDTH: 7.6 cm (3 in)

### PACKAGED DIMENSIONS

21.6 x 7.6 x 4.0 cm  
(8.5 x 3 x 1.6 in)

WT: 0.23 kg (0.5 lbs)